


August 12, 2010

NEW PRODUCT**M3PFHxS**

The global distribution of perfluoroalkylsulfonates has generated the need for analytical methods that facilitate the accurate analysis of these perfluorinated compounds in a variety of matrices. This class of environmental contaminants has been analyzed by both commercial and research laboratories for a number of years, but carbon-13 mass-labelled reference standards have only been available for perfluoro-1-octanesulfonate.

For this reason, **Wellington** has synthesized **M3PFHxS** in order to complement our existing selection of mass-labelled perfluoroalkylsulfonate reference standards.

**M3PFHxS****NEW**

Catalogue Number	Product (methanol)	Qty/Conc
MPFHxS	Sodium Perfluoro-1-hexane[¹⁸ O ₂]sulfonate	1.2 ml 50 µg/ml
M3PFHxS	Sodium Perfluoro-1-[1,2,3- ¹³ C ₃]hexanesulfonate	1.2 ml 50 µg/ml
MPFOS	Sodium Perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	1.2 ml 50 µg/ml
M8PFOS	Sodium Perfluoro-1-[¹³ C ₈]octanesulfonate	1.2 ml 50 µg/ml

Please contact your local distributor or info@well-labs.com for pricing and delivery.

Visit our website (www.well-labs.com) for a complete listing of our new products.

ISO 9001



CERT-0025901